

TUBE EVAC™

www.tubeevac.com

How to use TUBE EVAC™

1. Place the tubing in the groove of the tube-evac close to where it exits the skin.
2. With your other hand, hold the tubing firmly between the skin and the tube-evac.



3. Firmly squeeze the tube-evac and pull it towards the collection bulb.

4. Advance other hand as necessary to hold tubing.

5. The tube-evac rolls down the tubing very easily with a small but firm squeeze to evacuate small fibrinous material and or clots to maintain a functioning drain.



Lighthouse for the Blind
(800) 542-3697

TUBE-EVAC

Tube-Evac: A God inspired device created by the husband of a breast cancer patient. The device is designed to "clean" surgical tubing attached to a drain placed in a surgical site. It has no competition in the medical field.

Tube-Evac: Manufactured in USA under strict quality controls.

Tube-Evac: Patents: Design & Utility Patents Granted
International Patents both Pending & Granted

Tube-Evac: Designed to clear 7-10 & 15-19 French gauge tubing.

Tube-Evac: Prevents over-stretching and distortion of tubing.

Tube-Evac: Offers a quick, easy and safe way for doctors, nurses and patients to clear surgical tubing. Thousands have been used at major health care facilities.

Tube-Evac: Single patient use only. To be discarded once drains are discontinued.

Tube-Evac: Creates peace of mind for patient and caregiver with safe, effective ease of use.

Tube-Evac: Web site with demonstration video, manufacturer and order information available at tubeevac.com or Email: Customerservice@lhbindustries.com, Phone 1-800-542-3697, Fax 314-423-6918
Lighthouse for the Blind
10440 Trenton Avenue
St. Louis, MO 63132



**BOOTH 1663 at Plastic Surgery the Meeting
Stop by for a complimentary Tube-Evac device**